



InviLisa®

THE GOLD STANDARD IN ALLERGEN DETECTION: **INVILISA® ELISA KITS** FOR FOOD SAFETY

The InviLisa® ELISA kits provide a comprehensive solution for the quantification of food allergens from a variety of sources, including different fresh and processed food products, rinse waters and environmental swab samples.

ELISA KITS **HAVE NUMEROUS** **ADVANTAGES**

- Fast 60-minute (20:20:20) sandwich assay format
- Specific and highly sensitive detection of allergenic proteins
- Reliable quantitation for quality control
- Validated according to the AOAC guidelines for food allergens
- 96-well format with detachable strips (8x12)

FOR A COMPLETE **ALLERGEN TESTING** **SOLUTION**

Explore our PCR and LFA product range




InviScreen®



InviRapid®

InviScreen® Real-Time PCR Detection
InviRapid® Lateral Flow Test Strips



GET READY TO **SIMPLIFY YOUR WORKFLOW** WITH INVITEK DIAGNOSTICS



GET READY TO **SIMPLIFY YOUR WORKFLOW**
WITH INVITEK DIAGNOSTICS

Product	LOD	Quantitation Range	Catalogue Number
Inviliisa® Casein ELISA Kit	0.05 mg/kg	0.2 - 5 mg/kg	6032001200
Inviliisa® Milk ELISA Kit	0.1 mg/kg	0.5 - 10 mg/kg	6032002200
Inviliisa® Gluten ELISA Kit	0.1 mg/kg	2.5 - 50 mg/kg	6032003200
Inviliisa® Egg ELISA Kit	0.1 mg/kg	1.0 - 10 mg/kg	6032004200
Inviliisa® β-Lactoglobulin ELISA Kit	2.0 mg/kg	10 - 300 mg/kg	6032005200
Inviliisa® Soya ELISA Kit	0.3 mg/kg	0.8 - 10 mg/kg	6032006200

Standard package size 96 tests

CONTACT US

info@invitek.com

INVITEK DIAGNOSTICS PORTUGAL


ALS Life Sciences, SA
Zona Industrial de Tondela, ZIM II, Lote 6
3460-070 Tondela, Portugal

+351 232 817 817

INVITEK DIAGNOSTICS GERMANY

Invitek Molecular GmbH
Robert-Roessle-Str. 10
D-13125 Berlin, Germany

BE CURIOUS
AND SEE WHAT'S IN FOR YOU



invitek.com